

Product Datasheet

Human HSPb3 ELISA Kit (Ready to Use) (orb550995)

Description	Human Heat Shock Protein Beta 3 (HSPb3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q12988
Sensitivity	0.29ng/mL
Entrez	8988
Expiration Date	Please enquire.