

Product Datasheet

Human IHH ELISA Kit (Ready to Use) (orb551002)

Description	Human Hedgehog Homolog, Indian (IHH) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q14623
Sensitivity	0.057ng/mL
Entrez	3549
Expiration Date	Please enquire.