

Product Datasheet

Human Hemoglobin Beta ELISA Kit (Ready to Use) (orb551008)

Description	Human Hemoglobin Beta (Hemoglobin Beta) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Human
Range	1.56-100µg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at - 20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or erythrocyte lysates.
Sample Types	serum, plasma and erythrocyte lysates.
Uniprot ID	P68871
Sensitivity	0.45µg/mL
Entrez	3043
Expiration Date	Please enquire.