

Product Datasheet

Human HBg1 ELISA Kit (Ready to Use) (orb551009)

Description	Human Hemoglobin Gamma 1 (HBg1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for c...
Reactivity	Human
Range	468-30000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or erythrocyte lysates.
Sample Types	serum, plasma and erythrocyte lysates.
Uniprot ID	P69891
Sensitivity	122.6ng/mL
Entrez	3047
Expiration Date	Please enquire.