

Product Datasheet

Human HIF1a ELISA Kit (Ready to Use) (orb551088)

Description	Human Hypoxia Inducible Factor 1 Alpha (HIF1a) Sandwich ELISA Kit Ready to Use: Pre-diluted detection...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q16665
Sensitivity	0.055ng/mL
Entrez	3091
Expiration Date	Please enquire.