

Product Datasheet

Human IREM1 ELISA Kit (Ready to Use) (orb551093)

Description	Human Immune Receptor Expressed On Myeloid Cells 1 (IREM1) Sandwich ELISA KitReady to Use: Pre-dilut...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8TDQ1
Sensitivity	0.062ng/mL
Entrez	146722
Expiration Date	Please enquire.