

Product Datasheet

Human IGF2BP3 ELISA Kit (Ready to Use) (orb551118)

Description	Human Insulin Like Growth Factor 2 mRNA Binding Protein 3 (IGF2BP3) Sandwich ELISA KitReady to Use: ...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O00425
Sensitivity	0.114ng/mL
Entrez	10643
Expiration Date	Please enquire.