

Product Datasheet

Human IFI16 ELISA Kit (Ready to Use) (orb551165)

Description	Human Interferon Gamma Inducible Protein 16 (IFI16) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q16666
Sensitivity	0.049ng/mL
Entrez	3428
Expiration Date	Please enquire.