

Product Datasheet

Human IL1z ELISA Kit (Ready to Use) (orb551182)

Description	Human Interleukin 1 Zeta (IL1z) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for c...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9NZH6
Sensitivity	24pg/mL
Entrez	27178
Expiration Date	Please enquire.