

Product Datasheet

Human IL13Ra2 ELISA Kit (Ready to Use) (orb551192)

Description	Human Interleukin 13 Receptor Alpha 2 (IL13Ra2) Sandwich ELISA Kit Ready to Use: Pre-diluted detectio...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q14627
Sensitivity	0.129ng/mL
Entrez	3598
Expiration Date	Please enquire.