

Product Datasheet

Human IL17D ELISA Kit (Ready to Use) (orb551198)

Description	Human Interleukin 17D (IL17D) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for con...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8TAD2
Sensitivity	5.3pg/mL
Entrez	53342
Expiration Date	Please enquire.