

Product Datasheet

Human JKAMP ELISA Kit (Ready to Use) (orb551252)

Description	Human JNK1/MAPK8 Associated Membrane Protein (JKAMP) Sandwich ELISA KitReady to Use: Pre-diluted det...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Sensitivity	0.35ng/mL
Expiration Date	Please enquire.