

Product Datasheet

Human KLK10 ELISA Kit (Ready to Use) (orb551259)

Description	Human Kallikrein 10 (KLK10) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	043240
Sensitivity	0.048ng/mL
Entrez	5655
Expiration Date	Please enquire.