

Product Datasheet

Human KLK12 ELISA Kit (Ready to Use) (orb551261)

Description	Human Kallikrein 12 (KLK12) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Human
Range	15.6-1000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9UKR0
Sensitivity	7.14ng/mL
Entrez	43849
Expiration Date	Please enquire.