

Product Datasheet

Human KLK14 ELISA Kit (Ready to Use) (orb551263)

Description	Human Kallikrein 14 (KLK14) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9P0G3
Sensitivity	0.044ng/mL
Entrez	43847
Expiration Date	Please enquire.