

Product Datasheet

Human KLK5 ELISA Kit (Ready to Use) (orb551267)

Description	Human Kallikrein 5 (KLK5) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, saliva, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, saliva, tissue homogenates and other biological fluids.
Uniprot ID	Q9Y337
Sensitivity	29pg/mL
Entrez	25818
Expiration Date	Please enquire.