

Product Datasheet

Human KLF5 ELISA Kit (Ready to Use) (orb551316)

Description	Human Kruppel Like Factor 5, Intestinal (KLF5) Sandwich ELISA Kit Ready to Use: Pre-diluted detection...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q13887
Sensitivity	0.059ng/mL
Entrez	688
Expiration Date	Please enquire.