

## Product Datasheet

# Human LACRT ELISA Kit (Ready to Use) (orb551321)

<b>Description</b>	Human Lacritin (LACRT) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient...
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilitysaliva, tears or other biological fluids.
<b>Sample Types</b>	saliva, tears and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9GZZ8</a>
<b>Sensitivity</b>	0.055ng/mL
<b>Entrez</b>	<a href="#">90070</a>
<b>Expiration Date</b>	Please enquire.