

Product Datasheet

Human LAMR1 ELISA Kit (Ready to Use) (orb551343)

Description	Human Laminin Receptor 1 (LAMR1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for ...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P08865
Sensitivity	29pg/mL
Entrez	3921
Expiration Date	Please enquire.