

Product Datasheet

Human LCN1 ELISA Kit (Ready to Use) (orb551390)

Description	Human Lipocalin 1 (LCN1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, saliva, tears or other biological fluids.
Sample Types	serum, plasma, saliva, tears and other biological fluids.
Uniprot ID	P31025
Sensitivity	29pg/mL
Entrez	3933
Expiration Date	Please enquire.