

Product Datasheet

Human LRP6 ELISA Kit (Ready to Use) (orb551410)

Description	Human Low Density Lipoprotein Receptor Related Protein 6 (LRP6) Sandwich ELISA KitReady to Use: Pre-...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	075581
Sensitivity	0.059ng/mL
Entrez	4040
Expiration Date	Please enquire.