

Product Datasheet

Human MASP1 ELISA Kit (Ready to Use) (orb551460)

Description	Human Mannose Associated Serine Protease 1 (MASP1) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Human
Range	93.7-6000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P48740
Sensitivity	32.12pg/mL
Entrez	5648
Expiration Date	Please enquire.