
Product Datasheet

**IKK alpha/IKK beta (phospho-Ser180/Ser181)
antibody (orb5515)**

Description

IKK alpha/IKK beta (phospho-

Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IKK alpha/IKK beta around the phosphorylation site of Ser180/Ser181 (CT(p-S)FV)
Target	Phospho-IKK alpha (Ser180) + IKK beta (Ser181)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85/87 kDa
Uniprot ID	O15111
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.