

Product Datasheet

Human Melatonin Receptor 1A ELISA Kit (Ready to Use) (orb551503)

Description	Human Melatonin Receptor 1A (Melatonin Receptor 1A) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P48039
Sensitivity	0.219ng/mL
Entrez	4543
Expiration Date	Please enquire.