

## Product Datasheet

# Human MAP1LC3b ELISA Kit (Ready to Use) (orb551540)

<b>Description</b>	Human Microtubule Associated Protein 1 Light Chain 3 Beta (MAP1LC3b) Sandwich ELISA Kit Ready to Use:...
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Sensitivity</b>	0.068ng/mL
<b>Expiration Date</b>	Please enquire.