

## Product Datasheet

### Human MAOA ELISA Kit (Ready to Use) (orb551561)

<b>Description</b>	Human Monoamine Oxidase A (MAOA) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for ...
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Sensitivity</b>	0.114ng/mL
<b>Expiration Date</b>	Please enquire.