

Product Datasheet

Human MRP1 ELISA Kit (Ready to Use) (orb551574)

Description	Human Motility Related Protein (MRP1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents...
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P33527
Sensitivity	12.4pg/mL
Entrez	4363
Expiration Date	Please enquire.