

Product Datasheet

Human MAdCAM1 ELISA Kit (Ready to Use) (orb551589)

Description	Human Mucosal Addressin Cell Adhesion Molecule 1 (MAdCAM1) Sandwich ELISA KitReady to Use: Pre-dilut...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, urine or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, urine and other biological fluids.
Uniprot ID	Q13477
Sensitivity	0.125ng/mL
Entrez	8174
Expiration Date	Please enquire.