

Product Datasheet

Human MCL1 ELISA Kit (Ready to Use) (orb551602)

| | |
|--------------------------|---|
| Description | Human Myeloid Cell Leukemia Sequence 1, Bcl2 Related (MCL1) Sandwich ELISA KitReady to Use: Pre-dilu... |
| Reactivity | Human |
| Range | 0.312-20ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids. |
| Sample Types | tissue homogenates, cell lysates and other biological fluids. |
| Uniprot ID | Q07820 |
| Sensitivity | 0.118ng/mL |
| Entrez | 4170 |
| Expiration Date | Please enquire. |