

Product Datasheet

Human MYH3 ELISA Kit (Ready to Use) (orb551618)

Description	Human Myosin Heavy Chain 3, Skeletal Muscle, Embryonic (MYH3) Sandwich ELISA KitReady to Use: Pre-di...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P11055
Sensitivity	0.126ng/mL
Entrez	4621
Expiration Date	Please enquire.