

Product Datasheet

Human NAT8L ELISA Kit (Ready to Use) (orb551643)

Description	Human N-Acetyltransferase 8 Like Protein (NAT8L) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8N9F0
Sensitivity	0.064ng/mL
Entrez	339983
Expiration Date	Please enquire.