

Product Datasheet

Human NT-ProBNP ELISA Kit (Ready to Use) (orb551653)

Description	Human N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) Sandwich ELISA Kit Ready to Use: Pre-dilut...
Reactivity	Human
Range	39-2500pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	15pg/mL
Expiration Date	Please enquire.