

Product Datasheet

Human NT-ProCNP ELISA Kit (Ready to Use) (orb551654)

Description	Human N-Terminal Pro C-Type Natriuretic Peptide (NT-ProCNP) Sandwich ELISA KitReady to Use: Pre-dilu...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	4.42pg/mL
Expiration Date	Please enquire.