

Product Datasheet

Human NF2 ELISA Kit (Ready to Use) (orb551696)

| | |
|--------------------------|---|
| Description | Human Neurofibromin 2 (NF2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve... |
| Reactivity | Human |
| Range | 0.312-20ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids. |
| Sample Types | tissue homogenates and other biological fluids. |
| Uniprot ID | P35240 |
| Sensitivity | 0.128ng/mL |
| Entrez | 4771 |
| Expiration Date | Please enquire. |