

Product Datasheet

Human NEFH ELISA Kit (Ready to Use) (orb551698)

Description	Human Neurofilament, Heavy Polypeptide (NEFH) Sandwich ELISA Kit Ready to Use: Pre-diluted detection ...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P12036
Sensitivity	5.2pg/mL
Entrez	4744
Expiration Date	Please enquire.