

## Product Datasheet

### Human NEFL ELISA Kit (Ready to Use) (orb551699)

<b>Description</b>	Human Neurofilament, Light Polypeptide (NEFL) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
<b>Reactivity</b>	Human
<b>Range</b>	15.6-1000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, cerebrospinal fluid or other biological fluids.
<b>Sample Types</b>	serum, plasma, cerebrospinal fluid and other biological fluids.
<b>Uniprot ID</b>	<b>P07196</b>
<b>Sensitivity</b>	5.9pg/mL
<b>Entrez</b>	<b>4747</b>
<b>Expiration Date</b>	Please enquire.