

Product Datasheet

Human NT ELISA Kit (Ready to Use) (orb551724)

| | |
|--------------------------|---|
| Description | Human Neurotensin (NT) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient... |
| Reactivity | Human |
| Range | 4.68-300pg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | P30990 |
| Sensitivity | 1.35pg/mL |
| Entrez | 4922 |
| Expiration Date | Please enquire. |