

Product Datasheet

Human NT4 ELISA Kit (Ready to Use) (orb551728)

Description	Human Neurotrophin 4 (NT4) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P34130
Sensitivity	12.3pg/mL
Entrez	4909
Expiration Date	Please enquire.