

Product Datasheet

Human NEXN ELISA Kit (Ready to Use) (orb551734)

Description	Human Nexilin (NEXN) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient a...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q0ZGT2
Sensitivity	6.3pg/mL
Entrez	91624
Expiration Date	Please enquire.