

Product Datasheet

Human NOTCH2 ELISA Kit (Ready to Use) (orb551757)

Description	Human Notch Homolog 2 (NOTCH2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q04721
Sensitivity	0.33ng/mL
Entrez	4853
Expiration Date	Please enquire.