

Product Datasheet

Human NOTCH2NL ELISA Kit (Ready to Use) (orb551758)

Description	Human Notch2 N-Terminal Like Protein (NOTCH2NL) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q7Z3S9
Sensitivity	0.117ng/mL
Entrez	388677
Expiration Date	Please enquire.