

Product Datasheet

Human OSCAR ELISA Kit (Ready to Use) (orb551799)

Description	Human Osteoclast Associated Receptor (OSCAR) Sandwich ELISA Kit Ready to Use: Pre-diluted detection r...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8IYS5
Sensitivity	0.136ng/mL
Entrez	126014
Expiration Date	Please enquire.