

Product Datasheet

Human PP ELISA Kit (Ready to Use) (orb551824)

Description	Human Pancreatic Polypeptide (PP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for...
Reactivity	Human
Range	46.8-3000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01298
Sensitivity	13.3pg/mL
Entrez	5539
Expiration Date	Please enquire.