

Product Datasheet

Human PTH ELISA Kit (Ready to Use) (orb551833)

Description	Human Parathyroid Hormone (PTH) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for c...
Reactivity	Human
Range	1.56-100pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01270
Sensitivity	0.49pg/mL
Entrez	5741
Expiration Date	Please enquire.