

Product Datasheet

Human PMP22 ELISA Kit (Ready to Use) (orb551866)

Description	Human Peripheral Myelin Protein 22 (PMP22) Sandwich ELISA Kit Ready to Use: Pre-diluted detection rea...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q01453
Sensitivity	0.064ng/mL
Entrez	5376
Expiration Date	Please enquire.