

Product Datasheet

Human PXDN ELISA Kit (Ready to Use) (orb551870)

Description	Human Peroxidasin Homolog (PXDN) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for ...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q92626
Sensitivity	0.059ng/mL
Entrez	7837
Expiration Date	Please enquire.