

Product Datasheet

Human PRDX1 ELISA Kit (Ready to Use) (orb551871)

Description	Human Peroxiredoxin 1 (PRDX1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for con...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q06830
Sensitivity	0.126ng/mL
Entrez	5052
Expiration Date	Please enquire.