

## Product Datasheet

### Human PRDX5 ELISA Kit (Ready to Use) (orb551874)

<b>Description</b>	Human Peroxiredoxin 5 (PRDX5) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for con...
<b>Reactivity</b>	Human
<b>Range</b>	0.781-50ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability tissue homogenates, cerebrospinal fluid or other biological fluids.
<b>Sample Types</b>	tissue homogenates, cerebrospinal fluid and other biological fluids.
<b>Uniprot ID</b>	<a href="#">P30044</a>
<b>Sensitivity</b>	0.34ng/mL
<b>Entrez</b>	<a href="#">25824</a>
<b>Expiration Date</b>	Please enquire.