

Product Datasheet

Human PIK3C2a ELISA Kit (Ready to Use) (orb551909)

Description	Human Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide (PIK3C2a) Sandwich ELISA KitReady to Use: ...
Reactivity	Human
Range	3.12-200U/L
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O00443
Sensitivity	1.6U/L
Entrez	5286
Expiration Date	Please enquire.