

Product Datasheet

Human PDGFAB ELISA Kit (Ready to Use) (orb551953)

Description	Human Platelet Derived Growth Factor AB (PDGFAB) Sandwich ELISA Kit Ready to Use: Pre-diluted detecti...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	27pg/mL
Expiration Date	Please enquire.